

**MANUFACTURE PRODUCTS
WITH LOVE**

**SERVE CUSTOMERS
WITH HEART**

TIANJIN HUAHONG TECHNOLOGY CO., LTD.

Add: A01, Plant B, No.278, Hang Kong Road, Tianjin Pilot Free Trade Zone

(Air Port Industrial Park), Tianjin 300308, China

Tel: 86-22-58775333

Fax: 86-22-58775353

<http://www.hh-technology.com>

**Twist Lancet
Safety Lancet
Lancing Device
Insulin Pen Needle**

COMPANY PROFILE

Huahong has been developing, manufacturing, and marketing capillary blood sampling devices for over 10 years since 2010. We are the first enterprise in China who automates lancet production. Now, as one of the leading manufacturers of capillary blood sampling devices, Huahong's product has been distributed in many countries including USA, Canada, Europe, Korea, Japan, Brazil, Africa, etc.

With the mission of "people-oriented, create value", Huahong, continuously, innovates products as well as production processes and had obtained over 36 innovation patents and utility model patents.

Huahong is dedicated to providing health care workers and patients around the world with safe and convenient capillary blood sampling devices.

QUALITY FIRST

MDD 93/42/EEC

FDA Registered

EN ISO 13485: 2016

UL 21 CFR Part 820

CGMP & KGMP



Huahong Technology founded.
 Safety lancet and lancing device launched.
 Automate lancet injection.

2010

Safety lancet patented in China.
2012

Attained SA 8000 certificate.
2013

Publicly listed.
 Passed world-class pharmacy stores' third-party social responsibility audits.
2015

Safety lancet achieved award in Tianjin city government.
2016

The smallest safety lancet in the world launched
2017

Insulin pen needle launched.
2020

Safety lancet patented in USA.
2018

Lancing Device

- Unique patented design to minimize the pain level.
- Non-touch removal to control the risk of accidental injury.
- Compatible with most lancets in market.



UniShine® (Type HH-XVI-T)



Non-touch removal



UniShine® (Type HH-XVI-T)



Simple and elegant appearance design



UniVision® (Type HH-XIII-T)



10 depths available for different skin types




UniDolphin® (Type HH-X-T)



Window depth indicator


Elastic design



Invention Patent : ZL 2014 1 0314138.9

Interference between lancet and sleeve effectively reduces vibration of the needle, and makes needle go straight, which will help get maximum comfort.

Mini Lancing Device



ExtremeStar™ (Type HH-XIX)

Mini Design

Mini dimension, easy to use and carry.

Safety

Automatic retraction after use.

Convenience

Compatible with most lancet.

Various Depths

5 depths available for different skin types.

Code	Specification	Description	Color	Package
I16B015002	UniShine	10 depths	White/Black	1 pcs/bag 500 pcs/carton (or per customers' requests)
I13B015001	UniVision	5 depths	Gray+White	
I10B015001	UniDolphin	5 depths	Gray+White	
I10B015002	UniDolphin	10 depths	Gray+White	
I19B015001	ExtremeStar	5 depths	Gray/Black	

Twist Lancet

- Irradiation sterilization.
- 100,000-grade cleanroom.
- 21G-33G complete specifications.
- Compatible with most lancing devices.

TWIST LANCET




Tri-beveled needle tip



Transparent color



Straight needle inside



Stable exposed needle length



Compatible with most lancing devices

Twist Lancet 33G

- Large-scale production.
- More comfortable sampling experience.

TWIST LANCET



Needle diameter is only 0.21mm



Compatible with most lancing devices



Better choice for pain-sensitive user

Code	Gauge	Color	Package
IR21100200	21G	Light Blue	
IR23100200	23G	Light Blue	
IR26100200	26G	Light Blue	100 pcs/box
IR28100200	28G	Transparent Blue	200 boxes/carton
IR30100200	30G	Transparent Blue	(or per customers' requests)
IR33100200	33G	Light Blue	

Safety Lancet (Type XVII)

- Pressure activated.
- Tri-beveled needle tip to reduce the pain level.
- Automatic retraction for single use only.
- Irradiation sterilization and sealed needle to ensure 5 years of shelf life.



RuidyLance®



US&CN patented



Ergonomic handle design



Invention Patent

US: 9.826.927











CN: ZL 2012 1 0080546.3

- Stable firing system
- Elastic inner sleeve
- No repeated puncture

RuidyLance®—more than 10 years of guardianship

Low pain level
 Safe and reliable
 Convenient to use



Code	Specification	Color	Blood	Package
I17A2122100	21G*2.2mm			100 pcs/box 20 boxes/carton (or per customers' requests)
I17A2322100	23G*2.2mm			
I17A2618100	26G*1.8mm			
I17A2818100	28G*1.8mm			
I17A3016100	30G*1.6mm			

Safety Lancet (Type XXIII)

- Pressure activated.
- Tri-beveled needle tip to reduce the pain level.
- Automatic retraction for single use only.
- Irradiation sterilization and sealed needle to ensure 5 years of shelf life.

Get your Concise MejrLance®

Cost-effective
Safe and reliable
Convenient to use



MejrLance®











 4 components

 Ergonomic handle design

 Anti-slip texture to increase friction

 Firing system to be patented by US



Code	Specification	Color	Blood	Package
I23A2118100	21G*1.8mm			100 pcs/box 50 boxes/carton (or per customers' requests)
I23A2318100	23G*1.8mm			
I23A2616100	26G*1.6mm			
I23A2814100	28G*1.4mm			
I23A3014100	30G*1.4mm			

Mini Safety Lancet (Type XXV)

- Pressure activated.
- Tri-beveled needle tip to reduce the pain level.
- Automatic retraction for single use only.
- Irradiation sterilization and sealed needle to ensure 5 years of shelf life.



Mini RuidyLance®



US&CN patented



Easy to carry and use with mini dimension



Invention Patent

US: 9.826.927











CN: ZL 2012 1 0080546.3

- Stable firing system
- Elastic inner sleeve
- No repeated puncture

Get your Mini RuidyLance®

Eco-friendly
 Low pain level
 Elegant appearance



Code	Specification	Color	Blood	Package
I25A2122100	21G*2.2mm			100 pcs/box 40 boxes/carton (or per customers' requests)
I25A2322100	23G*2.2mm			
I25A2618100	26G*1.8mm			
I25A2818100	28G*1.8mm			
I25A3016100	30G*1.6mm			

Mini Safety Lancet (Type XXII)

- Pressure activated.
- Tri-beveled needle tip to reduce the pain level.
- Automatic retraction for single use only.
- Irradiation sterilization and sealed needle to ensure 5 years of shelf life.



SilkLance®



US&CN patented



Easy to carry and use with mini dimension



Invention Patent

US: 9.826.927

CN: ZL 2012 1 0080546.3

- Stable firing system
- Elastic inner sleeve
- No repeated puncture



Get your elegant SilkLance®

- Eco-friendly
- Low pain level
- Good choice for homecare



Code	Specification	Color	Blood	Package
I22A2122100	21G*2.2mm			100 pcs/box 40 boxes/carton (or per customers' requests)
I22A2318100	23G*1.8mm			
I22A2616100	26G*1.6mm			
I22A2816100	28G*1.6mm			
I22A3014100	30G*1.4mm			

Adjustable Safety Lancet (Type XXIV)

- Pressure activated.
- Tri-beveled needle tip to reduce the pain level.
- Automatic retraction for single use only.
- Irradiation sterilization and sealed needle to ensure 5 years of shelf life.




MejorLance® Adjustable Safety Lancet



Code	Specification	Depths	Color	Package
I24A2118100	21G	2.2mm 1.8mm 1.4mm		100 pcs/box 50 boxes/carton (or per customers' requests)
I24A2318100	23G	2.2mm 1.8mm 1.4mm		
I24A2616100	26G	2.0mm 1.6mm 1.2mm		
I24A2814100	28G	1.8mm 1.4mm 1.0mm		
I24A3014100	30G	1.8mm 1.4mm 1.0mm		

Top Fire Safety Lancet (Type XXVI)

- Push button to activate the lancet, which makes sampling more controllable.
- Ergonomic shape and grip handle make the sampling more comfortable.



Code	Specification	Color	Blood	Package
I26A2122100	21G*2.2mm			100 pcs/box 40 boxes/carton (Or per the customers' requests)
I26A2322100	23G*2.2mm			
I26A2618100	26G*1.8mm			
I26A2818100	28G*1.8mm			
I26A3016100	30G*1.6mm			

Heel Stick Safety Lancet (Type HHZ-II)

- Full specifications meet different blood sample collection needs.
- Irradiation sterilization and sealed individual package for 5 years of shelf life.



Qbutton®

↓ Only one step to get blood



Push down to release



Locked after use

Use Instructions



1 Select the test site.
 Clean the site with an antiseptic swab and allow to air dry.
 Remove from blister pack.




2 Place the heel stick against the baby's heel and activate by pressing downward on the white latch.




3 Discard it into a biohazard sharps container.

Motion blade design



Puncture blade



Motion blade

- Only reach the superficial blood vessels
- Enough blood sample
- Low pain level
- Quick wound healing



Code	Specification	Color	Package
IH06505004	0.65 mm *1.5 mm		50 pcs/box
IH08505004	0.85mm*1.75 mm		20 boxes/carton
IH10005004	1.0mm*2.5mm		(or per customers' requests)
IH15005004	1.5mm*2.8mm		


HEEL STICK

HEEL STICK

Insulin Pen Needle (Ordinary Type)

Comfortable

- Tri-beveled and precision sharpened needle tips.
- Thin wall needle to meet ISO9626 for optimal insulin flow.
- Advanced siliconized needle lubrication to reduce pain.



Universal

- Compatible with most insulin pens, such as
Lilly: HumaPen Luxura, HumaPen Memoir, HumaPen Ergo II, KwikPen.
Nova Nordisk: NovoPen Echo, NovoPen-4, FlexPen.
Owen Mumford: AutoPen.
Sanofi Aventis: ClikStar, SoloStar.
Amylin: SymlinPen.
- Please call Huahong for other pens.



Safe

- Transparent appearance, the needle is clear to observe.
- With Plasma treatment, needle tube and needle hub are bonded solidly.



PrickLess®—Designed for Maximum Comfort



Code	Specification	Color	Package
IP300510050	30G *5mm		100 pcs/box 50 boxes/carton (or per the customers' request)
IP300610050	30G *6mm		
IP300810050	30G *8mm		
IP310510050	31G *5mm		
IP310610050	31G *6mm		
IP310810050	31G *8mm		
IP320410050	32G *4mm		
IP320510050	32G *5mm		
IP320610050	32G *6mm		
IP320810050	32G *8mm		

Insulin Pen Needle (Safety Type)

PrickLess[®]—Designed for Maximum Comfort

Safety

Red sign provide visual confirmation that the pen needle is activated.
 Shield automatically locks after injection to avoid accidental injury.

Comfort

Needle is hidden from view before and after use for a less intimidating experience.



Flexibility

Flexibility to work with most insulin pens on the market.



Red sign shows it is used

Cannot be loaded again

Specification	Color	Package
30G*6mm		100 pcs/box 30 boxes/carton (or per customers' requests)
30G*8mm		
31G*6mm		
31G*8mm		

Notes

INSULIN PEN NEEDLE