



F.L. MEDICAL s.r.l. Unipersonale
 Via Enrico Mattei, 20 – 35038 TORREGLIA (Padova) – Italy
 Tel. +39 049 5211475 / +39 049 5212835 - Fax +39 049 5212566
 e-mail: info@flmedical.com - web site: www.flmedical.com
 C.F. e P.IVA 01134840287 - Cap Soc. 90.000 € i.v.
 Reg. Imp. di Padova n. 21695 - R.E.A. di Padova n. 187254

FICHA TÉCNICA / TECHNICAL INFORMATION
CÓDIGO / ITEM BK- 21014

TUBO 12x75 mm 5 ml CILÍNDRICO, SIN BORDE, EN POLIPROPILENO, BOLSA DE 250 PIEZAS
 TEST TUBE 12x75 mm 5 ml CYLINDRICAL, WITHOUT RIM, IN POLYPROPYLENE SORWALL TYPE IN
 BAG OF 250 PCS

Fabricante / Manufacturer	F.L. MEDICAL s.r.l.
Dirección / Address	Via Enrico Mattei, 20 – 35038 TORREGLIA (Padova) – Italy

1. INFORMACIÓN TÉCNICA / TECHNICAL INFORMATION:

CÓDIGO DEL PRODUCTO / CODE:	BK-21014
CAPACIDAD / CAPACITY:	5 ml
MATERIAL TUBO / CUP MATERIAL:	Polipropileno / Polypropylene
MATERIAL TAPÓN / CAP MATERIAL:	Polietileno / Polyethylene
APARIENCIA / APPEARANCE:	Cilíndrico / Cylindrical
DIMENSIONES EXTERNAS / OUTSIDE SIZES:	\varnothing inferior / \varnothing lower: 12 mm \varnothing superior / \varnothing higher: 12 mm Altura sin tapa / Height without cap: 75 mm Altura con tapa / Height with cap: N.A
COLOR TUBO / CUP COLOUR:	Neutral, transparente / Neutral, transparent
COLOR TAPÓN / CAP COLOUR:	No applicable / Not applicable
OTRAS CARACTERÍSTICAS / OTHER CHARACTERISTICS:	No applicable / Not applicable Superficie de escritura para datos paciente / frosted label for patient data: NAME, N°, DATE

2. ESTERILIZACIÓN / STERILIZATION:	No applicable / Not applicable
3. CANTIDAD POR CAJA / BOX QUANTITY:	Caja de 4000 pz / 4000 pcs. per box
4. CANTIDAD INTERIOR / INNER QUANTITY:	Bolsa de 250 piezas / Bag of 250 pcs
5. VOLUMEN CAJA / BOX VOLUME:	mc / cbm: 0,086 Medidas caja / box sizes: 59 x 38.5 x 38 cm
6. PESO CAJA / BOX WEIGHT:	Peso bruto / Gross weight: 6,0 Kg
7. INFORMACIÓN ADICIONAL / ADDITIONAL INFORMATION:	Producto IVD, marcado CE / IVD item, CE marked
8. ALMACENAMIENTO / STORAGE:	No ponerlos en contacto o cerca de llamas, evítese su liberación al medio ambiente ya que NO son biodegradables / Do not expose or leave near open flames, do not leave the product around as it is not biodegradable.